



# Hidden CATLAB

Integrated Microreactor-Mass Spectrometer  
for Reaction Testing, TPD/R/O and  
Pulse Chemisorption

# CATLAB Modules




# Standard CATLAB Layout



Module 2  
Mass Spectrometer

Module 1  
Microreactor and TCU/GCU

# CATLAB Schematic



# CATLAB Options


- Corrosion resistant upgrade - includes upgrade of one gas flow line and corrosive gas resistant sample line seals, and a gas dilution/purge valve mass spectrometer pumping system.
- Additional corrosion resistant feed lines if required.
- Additional 4 channel gas control unit integrated to provide 8 stream gas selection in total.
- Option to choose the maximum flow of each MFCs.
- Vapour Generation (gives vapour pressure equivalent to max liquid temperature ~30°C). Other options available for higher flow rates.
- 300 amu mass range option. Mass range to 1000 amu available for specialised applications

# QIC Series Gas Analysers




- 200 amu mass range, dual Faraday/ Electron Multiplier detector.
- Detection capability from 100 % to 0.1ppm.
- Fast scan speeds of 100 amu/s.
- < 500 ms response time to changes in gas concentrations.
- Low dead volume, flexible heated inlet capillary for fast response to gases and vapours.
- Soft ionisation for analysis of complex mixtures and organics.
- Can be used as a stand-alone gas analysis system or in combination with other equipment i.e. Thermal Analysis.

# CATLAB Control Software



- Control of Mass Spectrometer, Temperature and Gas Flows in one software package



# Temperature Control



- Multi-stage temperature ramps




# Gas Mixing Control



- MFC Gas Mix Library Editor


- Gas Mix Calculator

# Gas Mixing Control



- Gas switching controlled by either temperature or time

# MS Analysis Control



- Optimised multistage analysis - configure different analysis for different parts of the experiment


# Mass Spectrometer Control

**MID Setup**

Gas 1 Argon	40	20	36	38	0	0	0	0	0	0
	999	107	3	1	0	0	0	0	0	0
Gas 2 Hydrogen	2	1	0	0	0	0	0	0	0	0
	999	100	0	0	0	0	0	0	0	0
Gas 3 Water	18	17	16	2	20	19	0	0	0	0
	999	230	11	7	3	1	0	0	0	0
Gas 4 Oxygen	32	16	34	33	0	0	0	0	0	0
	999	114	4	1	0	0	0	0	0	0
Gas 5	0	0	0	0	0	0	0	0	0	0
	0	0	0	0	0	0	0	0	0	0
Gas 6	0	0	0	0	0	0	0	0	0	0
	0	0	0	0	0	0	0	0	0	0
Gas 7	0	0	0	0	0	0	0	0	0	0
	0	0	0	0	0	0	0	0	0	0
Gas 8	0	0	0	0	0	0	0	0	0	0
	0	0	0	0	0	0	0	0	0	0


- MID Scan setup for known species
- Automatic overlap removal

# Mass Spectrometer Control




- Multiple Bar scans can be configured for optimised sampling of unknowns

# MS Data Display




- Data plotted with x-axis as time or temperature

# MS Data Display



- 3D Bar Graph mode for easy identification of bar mode trends


# Data Analysis



- Peak fitting analysis routines
- Integrated area
- Baseline subtraction




# Pulse Chemisorption



- Injection of single gases or multiple reactants

# Pulse Chemisorption Quantification



- Uptake measurements
- Dispersion calculation
- Metal surface area
- Pulse Adsorption Isotherm

# Summary

- Integrated microreactor and mass spectrometer
- Both microreactor and mass spectrometer manufactured by Hiden Analytical
- Single integrated software package to control MS and microreactor
- Ideal for catalyst characterisation and reaction testing
- Designed for optimum mass spectrometer performance



## Selected Academic References

- **Electronic and Geometric Structure of Ce<sup>3+</sup> Forming Under Reducing Conditions in Shaped Ceria Nanoparticles Promoted by Platinum** . *J. Phys. Chem. C*, 2014, 118 (4), pp 1974–1982. O. V. Safonova, A. Guda, C. Paun, N. Smolentsev, P. M. Abdala, G. Smolentsev, M. Nachtegaal, J. Szlachetko, M. A. Soldatov, A. V. Soldatov, and J. A. van Bokhoven
- **The effect of coke deposition on the activity and selectivity of HZSM-5 zeolite during ethylbenzene alkylation reaction in the presence of ethanol.** *Catal. Sci. Technol.*, 2014, **4**, 1017-1027. S. S. Khattaf, C. D'Agostino, M. N. Akhtar, N. O. Al-Yassir, N. Tan and L. F. Gladden.
- **Reaction Kinetics of C<sub>3</sub>H<sub>6</sub> Oxidation for Various Reaction Pathways Over Diesel Oxidation Catalysts** . *Topics in Catalysis*, December 2013, Volume 56, Issue 18-20, pp 1916-1921. H. Oh, I. S. Pieta, J. Luo, W. S. Epling
- **An investigation of the role of surface nitrate species in the oxidation of propene on a Pt-based diesel oxidation catalyst..** *Catal. Sci. Technol.*, 2013, **3**, 2349-2356 , S. Chansai, R. Burch, C. Hardacre, H. Oh, and W. S. Epling.
- **Molybdenum oxide on Fe<sub>2</sub>O<sub>3</sub> core-shell catalysts; probing the nature of the structural motifs responsible for methanol oxidation catalysts** . *ACS Catal.*, **2014**, 4 (1), pp 243–250, C. Brookes , P. P Wells , G. Cibin , N. Dimitratos , W. Jones , D. J. Morgan , and M. Bowker.
- **Chiral Co(II) metal-organic framework in the heterogeneous catalytic oxidation of alkenes under aerobic and anaerobic conditions.** *ACS Catal.*, 2014, 4, pp 1032–1039 Giulia Tuci, Giuliano Giambastiani, Stephanie Kwon, Peter Curran Stair, Randall Q. Snurr and Andrea Rossin

# Selected CATLAB Users

**PDVSA****UCL**  
Université  
catholique  
de Louvain**ETH** zürich

- Johnson Matthey
- Texas A&M University
- Hong Kong University
- Cambridge University
- Bayreuth University
- Research Complex at Harwell
- Sao Paolo University
- Bulgarian Academy of Sciences
- Université Catholique de Louvain
- Kunming University
- ETH Zurich
- PDVSA



Johnson Matthey



The University of Hong Kong

UNIVERSITY OF  
CAMBRIDGE

- 
- A photograph of a modern, two-story office building with a grey facade and large glass windows. The building has a prominent entrance on the left side. A large, semi-transparent white circle is overlaid on the left side of the image, containing text. The sky is clear blue, and there are some trees and bushes in the foreground.
- [www.HidenAnalytical.com](http://www.HidenAnalytical.com)
  - The Hiden website is an excellent resource with product pages, brochures, catalogues, product pages with some application notes, presentation and other information.
  - Contact +44 1925 445225 for direct support.