

Plunger Sampling Pipette acc. to Hensen

Easy retrieval of precise small plankton sub-samples

ORDERING INFORMATION	
Products	
435 090	Plunger Sampling Pipette acc. to Hensen , capacity 0.1 ml
435 091	Plunger Sampling Pipette acc. to Hensen , capacity 0.25 ml
435 092	Plunger Sampling Pipette acc. to Hensen , capacity 0.5 ml
435 093	Plunger Sampling Pipette acc. to Hensen , capacity 1.0 ml
435 094	Plunger Sampling Pipette acc. to Hensen , capacity 2.5 ml
435 095	Plunger Sampling Pipette acc. to Hensen , capacity 5.0 ml
Spare Parts	
435 098	Spare Dilution Bottle , 250 ml capacity, made of Duranglass 50 (package of 10 pcs.)

ABOUT

The Plunger Sampling Pipettes acc. to Hensen are simple but precise tools to retrieve small plankton sub-samples. The volume of the sub-samples is defined by the accurately machined geometry of the piston, avoiding any volume errors. Each pipette is delivered completely with dilution bottle and stopper.



Folsom's Plankton Sample Divider

Easy division of large plankton samples



ABOUT

Folsom's Plankton Sample Divider provides a fast and easy method to divide large plankton samples into an amount suitable for examination. This instrument divides samples in halves. Simply repeat the operation until the desired amount is left. The transparent plastic drum of the divider is suitable for samples up to 100ml. The base is provided with a water level and two levelling screws for easy horizontal alignment.



ORDERING INFORMATION	
Products	
435 100	Folsom's Plankton Sample Divider