



HYDRO-BIOS

Smart Sampling

Datasheet Kolkwitz chamber 0.5 ml (435 015)



Plankton Chamber acc. to Kolkwitz 0.5 ml

Tray for easy enumeration of small amounts of phytoplankton.

The Plankton Chamber acc. to Kolkwitz has been designed for easy enumeration of small amounts of phytoplankton. Due to the small dimensions the chambers can be used with standard microscopes, inverted microscopes and stereo microscopes. Made from glass the chambers are precisely adjusted to their nominal volume. The cover plates are marked with an orthogonal counting grid (1 mm) for an easy zoning of the sample.

TECHNICAL DETAILS

length	33 mm
width	33 mm
height	1.3 mm
sample volume	0.5 ml
materials	glass
product ID	435 015



CHARACTERISTICS

- acc. to Kolkwitz