



HYDRO-BIOS

Smart Sampling

Datasheet chamber 12 x 2.4 ml (435 011 B)



counting chamber acc. to Gollasch 12 x 2.4 ml

Based on the Bogorov design

The Counting Chambers for Zooplankton have been designed for easy enumeration of large amounts of zooplankton. A plexiglass block with precisely machined grooves and a polished bottom provides best possible transparency. Due to the small dimensions all chambers can be used with stereo microscopes. The Ballast Water models, modified by Gollasch, are manufactured with grooves in line pattern, bisect line pattern and grid pattern. This very much alleviates the use of the chambers on moving vessels even under choppy sea conditions.

DESCRIPTION

+++ do not use alcohol for cleaning acrylic parts +++

TECHNICAL DETAILS

length	80 mm
width	100 mm
sample volume	12 x 2.4 ml
product ID	435 011 B



CHARACTERISTICS

- acc. to Gollasch