



HYDRO-BIOS

Smart Sampling

Datasheet chamber 40 x 70 mm (435 010)



Counting Chamber 40 x 70 mm

Based on the Bogorov design

The Counting Chambers for Zooplankton have been designed for easy enumeration of large amounts of zooplankton. A plexiglass block with precisely machined grooves and a polished bottom provides best possible transparency. Due to the small dimensions all chambers can be used with stereo microscopes. The standard models are manufactured with grooves in meander profile.

DESCRIPTION

+++ do not use alcohol for cleaning acrylic parts +++

TECHNICAL DETAILS

length	40 mm
width	70 mm
sample volume	9 ml
product ID	435 010



CHARACTERISTICS

- Based on Bogorov design