



HYDRO-BIOS

Smart Sampling

Datasheet chamber 80 x 100 mm (435 011)



Counting Chamber 80 x 100 mm

Based on the Bogorov design

The Counting Chambers for Zooplankton have been designed for easy enumeration of large amounts of zooplankton. A plexiglass block with precisely machined grooves and a polished bottom provides best possible transparency. Due to the small dimensions all chambers can be used with stereo microscopes. The standard models are manufactured with grooves in meander profile.

DESCRIPTION

+++ do not use alcohol for cleaning acrylic parts +++

TECHNICAL DETAILS

| | |
|----------------------|---------|
| length | 80 mm |
| width | 100 mm |
| sample volume | 22 ml |
| product ID | 435 011 |



CHARACTERISTICS

- Based on Bogorov design