



HYDRO-BIOS

Smart Sampling

Datasheet chamber 80 x 140 mm (435 012)



Counting Chamber 80 x 140 mm, 70ml

Based on the Bogorov design

The Counting Chambers for Zooplankton have been designed for easy enumeration of large amounts of zooplankton. A plexiglass block with precisely machined grooves and a polished bottom provides best possible transparency. Due to the small dimensions all chambers can be used with stereo microscopes. The standard models are manufactured with grooves in meander profile.

DESCRIPTION

+++ do not use alcohol for cleaning acrylic parts +++

TECHNICAL DETAILS

length	80 mm
width	140 mm
sample volume	70 ml
product ID	435 012



CHARACTERISTICS

- Based on Bogorov design