



HYDRO-BIOS

Smart Sampling

Datasheet Hensen plunger 0.1 ml (435 090)

Plunger Sampling Pipette acc. to Hensen 0.1 ml

Easy retrieval of precise small plankton sub-samples

The Plunger Sampling Pipettes acc. to Hensen are simple but precise tools to retrieve small plankton sub-samples. The volume of the sub-samples is defined by the accurately machined geometry of the piston, avoiding any volume errors. Each pipette is delivered completely with dilution bottle and stopper.



TECHNICAL DETAILS

sample volume	0.1 ml
materials	Glas, POM, stainless steel
product ID	435 090