



# HYDRO-BIOS

Smart Sampling

Datasheet Hensen plunger 0.5 ml ( 435 092 )

## Plunger Sampling Pipette acc. to Hensen 0.5 ml

Easy retrieval of precise small plankton sub-samples.

The Plunger Sampling Pipettes acc. to Hensen are simple but precise tools to retrieve small plankton sub-samples. The volume of the sub-samples is defined by the accurately machined geometry of the piston, avoiding any volume errors. Each pipette is delivered completely with dilution bottle and stopper.



## TECHNICAL DETAILS

<b>sample volume</b>	0.5 ml
<b>materials</b>	Glas, POM, stainless steel
<b>product ID</b>	435 092