




HYDRO-BIOS

Smart Sampling

Datasheet V-Fin Depressor 45kg V2A (438 419)

V-Fin Depressor 45kg V2A

The V-Fin Depressor gives negative lift to the towing rope during horizontal or oblique plankton net operation. This also stabilises the net in the water and allows for good sampling results. It consists of a lead-weighted body and a stainless steel rope with shackle.



TECHNICAL DETAILS

length	620mm
width	670mm
height	400mm
weight on air	45kg
materials	stainless steel, lead
product ID	438 419