



HYDRO-BIOS

Smart Sampling

Datasheet Rod closing mech. (437 320)

Rod-Held closing mechanism

Special Closing Mechanism

for Bottom Samplers acc. to Ekman-Birge and Lenz. In shallow waters the Bottom Samplers acc. to Ekman-Birge and Lenz can easily be operated with the Rod-Held Closing Mechanism. The samplers will be pushed by hand into the water bed and then released by the closing mechanism of the sampler. Any rod with a diameter of 20 mm can be used.



TECHNICAL DETAILS

weight on air	0.3 kg
product ID	437 320